




Thin Film XRF Reference Samples RF4-100-S1749 and RF4-200-S1749



Sample dimensions:

The reference samples are available in two different designs: Fixed to a circular PEEK holder for easier handling or on a small silicon frame. The frame size is 10x10 mm² with a 100 nm or 200 nm thin usable area of 5x5 mm² in the center. The plastic holder is 3 mm thick with an outer diameter of 30 mm or 49 mm.

Element content:


The reference samples contain the 7 elements **Pb, La, Pd, Mo, Cu, Fe** and **Ca** plus **Si** and **N** from the silicon nitride membrane. Further, they may contain traces of C, Ar or other contents of ambient air that are not important for most XRF measurements.

Possible applications:

- Absorption free standard: no matrix correction necessary
- Mass depositions in the range of ng/mm² (1-3 atomic layers) permit quantification without the need to interpolate from higher values
- Wide selection of non-overlapping XRF lines, exact calibration curve with many points over a large energy range
- Adjustment of confocal μ -XRF
- Transmission experiments due to low thicknesses of 100 nm or 200 nm



Thin Film XRF Reference Samples RF4-100-S1749 and RF4-200-S1749



Energy spectrum of the 7-element reference sample RF4-200-S1749 measured at 40 keV excitation. The energy range from ~2 keV to ~40 keV is covered with peaks of comparable intensity.

Mass deposition:



Mass depositions on the sample are in the range from 1 ng/mm² to 20 ng/mm² and listed below for all elements contained.

The mass deposition values listed here are average values measured by AAS, ICP-OES and μ -XRF. Despite the very precise measurements these reference samples are no "Certified Reference Materials (CRM)".

Element	Mass (ng/mm ²)	Emission Lines (eV)	
		K α	L α
Pb	7.61 \pm 0.96	85335	10541
La	11.01 \pm 0.62	33298	4649
Pd	1.8 \pm 1.0	21123	2838
Mo	1.32 \pm 0.40	17444	2293
Cu	2.84 \pm 0.35	8040	930
Fe	5.04 \pm 0.87	6401	747
Ca	19.31 \pm 1.10	3691	341
Si	Substrate	1740	

Lateral homogeneity:

The lateral homogeneity of all elements deposited on these reference samples has been tested with μ -XRF mappings. The deviation is smaller than 1% over the entire sample area.

	Large area map SF1	μ beam "mapping" S10
Energy	26 keV	9.5 keV
Area	15 x 15 mm ²	1.2 x 1.2 mm ²
Beam size	0.8 x 0.4 mm ²	2.8 x 12 μ m ²
Step size	0.8 x 0.4 mm ²	~ 30 x 30 μ m ²
Cu K α		
La L α	